12ml Micro-Berty

Catalytic Packless Reactor:

Volume: 12 ml

Vessel MAWP*: 5,000 psi @ 650°F (345 bar @ 343°C)

5,000 psi @ 1,000°F (345 bar @ 538°C)

Material of Construction: Hastelloy® C276 / Body, Housing, Basket

* Maximum Allowable Working Pressure



Principle of Operation:


The Parker Autoclave Engineers' 12 ml Catalytic Packless Reactor with Berty internals is designed to provide all the benefits of small scale chemical research. It allows the research scientist to work with small quantities of catalyst and feedstocks which may be expensive and/or limited in availability. Reduced volumes are safer to work with and minimize waste disposal. In addition, the reduced volume of the 12 ml capacity has reduced dead volume, which minimizes side thermal reactions. The role of the catalyst thus may be more clearly tested.


The catalyst basket design lends itself to easy loading and discharge of catalyst pellets. It is patterned after the larger proven Berty basket reactor currently available in 3" and 5" capacity vessels. Facilitating tomorrows requirements in a small proven stirred reactor package for University and Research facilities the world over.

Features:

- Maximum blower speed of 5,000 rpm
- Operating pressures as high as 5,000 psi @ 1,000°F (345 bar @ 538°C)
- Hydrostatic test pressures: 7,000 psi @ room temperature, 8,700 psi @ room temperature (Hi-Temp Units)

Dimensions:









Standard Equipment:

The micro catalytic reactor internal parts are designed to be used in a standard MicroClave[®] assembly with minimal changes. A complete assembly is available as a standard unit or a retrofit package can be purchased to convert an existing MicroClave[®] into a catalytic reactor.

Micro Reactor Vessel: Confined gasket closure employs flange nut to lock body and housing.

Body, housing, and flange nut are Hastelloy C-276

Sealing Gasket: Confined gasket of silver-plated Inconel-X is designed for temperatures to 650°F.

The high temperature unit uses a gold-plated metal seal rated at 1,000°F

Cover: Cover is integral with the MagneDrive housing.

Capacity: Usable capacity with reactor internals in place is 12ml

Connection Collar Standard openings include:

Openings: • One 1/8" connection for safety head and pressure gauge

One thermowellOne sample tube

• Two inlet/outlet charging connections

Body Openings: • Two bottom connections: Thermocouples

Purge Connection: 1/8" (SW125) gas connection at top of MagneDrive allows for introduction of gas into

the vessel.

Pressure Gauge: Constant reading gauge has 2-1/2" diameter dial with Monel Bourdon tube.

Dual face dial reads 0-7500 psi and 517 bar.

Safety Head Assembly: Hastelloy-C safety head uses 3/16" flat rupture disc rated 4,781-5,000 psig @ 72°F, with

1/8" NPT female vent connection through top of bench stand to atmosphere.

Furnace: External band-type electric furnace. 120 or 240VAC single phase.

Cooling Coil: External cooling coil can be used for water or air cooling. Type 304 SS coil permits rapid

vessel cooling and temperature control.

MagneDrive Packless The AE Micro-Berty features a packless MagneDrive system. Rare earth magnetics provide

Drive System: high torque mixing capability. Packless magnetic-drive system eliminates leakage,

contamination and packing heat generation problems of conventional mixers. It provides continuous high speed rotary agitation without the danger of leakage or the downtime to

change worn packing.


Mixing System: MagneDrive rotary impeller system. Static torque 6 in lbs; net mixing horsepower 0.5 @

5,000 rpm. Special impeller for maximum dispersion.

Available with electric Electric: Variable speed rated 1/4 HP @ 3,450 rpm. Type TENV GP with 120V or 240V.

or air motor: Air: HP to 0.5 @ 6,000 rpm. Required air pressure 80 psig @ 27 cfm maximum.

12ml Catalytic Packless Reactor with Micro-Berty Basket:



Technical Specifications:

This internal recycle reactor is designed with a fixed, circular, screened catalyst bed and a top mounted vane type blower. Fluid circulation is directed downward along the vessel wall and deflected upward through the catalyst bed.

Reactants: Gas/Solids. Vapor/Solids

Typical Reactions: Oxidation, hydrogenation, catalyst testing

Basket Screen: 50x50 mesh, 0.009" (0.23 mm) wire and a nominal opening size of 0.011" (0.28 mm)

Inside Diameter: 1" (25.4 mm)

Basket Volume: 0.22 in.3 (3.6 cm³)

Free Volume: 0.94 in.3 (15.4 cm³)

Maximum Allowable 5,000 psig (345 bar)

Working Pressure:

Maximum Blower Speed: 5,000 RPM

Version: High Temperature - 1,000°F (538°C)

Catalog Number Prefix: CRBHT


Standard Temperature - 650°F (343°C) Catalog Number Prefix: CRB2HC

Common Customization: Special wire mesh size, special materials,

specific pressure/temperature ratings.
ASME code stamp (or CE mark for pressure

equipment Directive)

Standard Material: Hastelloy® C-276



Ordering Guide:

The following reactor assemblies include motor, thermocouples, and electrically heated 1,400°F (760°C) maximum furnace (for the voltage specified in the table). Be advised, motor controls, tachometer display, furnace controls and the display for the thermocouple are purchased as separate items. The Specifications and descriptions found in the drawings referenced in the table below supercede the specification information found in this guide. Consult factory for more information.

| Catalog Number | Description SS=ANSI 316 Stainless Steel HC=Hastelloy [®] C-276 | Motor | Power Source | Temperature Rating | General Arrangement Drawing Number | Reactor Subassembly Drawing Number | Weight lbs. |
|-------------------|---|-------|-----------------|-----------------------|--|--|-------------|
| CRB2HC05ZH16A | Micro Berty Reactor 12 cc HC | Air | 120V | 650°F (343°C) | 40A-3174 | 40A-3173 | 47 |
| CRB2HC05ZH16D | Micro Berty Reactor 12 cc HC | DC | 120V | 650°F (343°C) | 40A-3175 | 40A-3173 | 62 |
| CRBHT2HC05ZH16D | Micro Berty Reactor 12 cc HC | DC | 120V | 1000°F (538°C) | 40A-7791 | - | 62 |
| CRB2HC05ZH26A | Micro Berty Reactor 12 cc HC | Air | 240V | 650°F (343°C) | 40A-3174 | 40A-3173 | 47 |
| CRB2HC05ZH26D | Micro Berty Reactor 12 cc HC | DC | 240V | 650°F (343°C) | 40A-3175 | 40A-3173 | 62 |
| CRBHT2HC05ZH26D | Micro Berty Reactor 12 cc HC | DC | 240V | 1000°F (538°C) | 40A-7791 | - | 62 |

NOTES:

The circulating pressure generated by the impellers in the "Micro Series" reactors is low.

Parker Autoclave Engineers makes no claims about the ability to scale-up or correlate "Micro Series" catalytic reactors with any other process equipment.

WARNING

FAILURE, IMPROPER SELECTION OR IMPROPER USE OF THE PRODUCTS AND/OR SYSTEMS DESCRIBED HEREIN OR RELATED ITEMS CAN CAUSE DEATH, PERSONAL INJURY AND PROPERTY DAMAGE. This document and other information from Parker Hannifin Corporation, its subsidiaries and authorized distributors provide product and/or system options for further investigation by users having technical expertise. It is important that you analyze all aspects of your application and review the information concerning the product or system in the current product catalog. Due to the variety of operating conditions and applications for these products or systems, the user, through its own analysis and testing, is solely responsible for making the final selection of the products and systems and assuring that all performance, safety and warning requirements of the application are met. The products described herein, including without limitation, product features, specifications, designs, availability and pricing, are subject to change by Parker Hannifin Corporation and its subsidiaries at any time without notice.

Offer of Sale

The items described in this document are available for sale by Parker Hannifin Corporation, its subsidiaries or its authorized distributors. Any sale contract entered by Parker will be governed by the provisions stated in Parker's standard terms and conditions of sale (copy available upon request).

©2016 Parker Hannifin Corporation | Autoclave Engineers is a registered trademark of the Parker Hannifin Corporation

06-1143SE April2016





Instrumentation Products Division

Autoclave Engineers Operation 8325 Hessinger Drive Erie, PA 16509-4679 Tel: 814 860 5700 • Fax: 814 860 5718 www.AutoclaveEngineers.com **Caution!** Do not mix or interchange component parts or tubing with those of other manufacturers. Doing so is unsafe and will void warranty.

Caution! Parker Autoclave Engineers Valves, Fittings, Systems, and Tools are not designed to interface with common commercial instrument tubing and are designed to only connect with tubing manufactured to Parker Autoclave Engineers AES specifications. Failure to do so is unsale and will void warranty.

Bulletin CR-MICROBERTY

ENGINEERING YOUR SUCCESS.