
O-ring Closure

Pressure Vessel

Volume: 100 ml to 2000 ml

Vessel MAWP*: 10,500 to 60,000 psi @ 72°F (724 to 4137 bar @ 22°C) Material of Construction: 316 Stainless Steel, Alloy Steel, A-286



Principle of Operation:

Parker Autoclave Engineers O-ring closure pressure vessels, model OR and O-ring closure compression cylinders, model OC and OD, offer a simple, reliable and economical alternative to the Parker Autoclave Engineers self-sealing closure where high temperature operation is not required. Standard O-ring material is Nitrile; however, optional O-ring materials are available on special order. Parker Autoclave Engineers welcomes your special product request. Please contact the factory for review of your application.

The main nut screws into the vessel body until it tightens the cover against the angle in the body, as shown in the figure. The groove in the cover is sized to put the proper squeeze on the O-ring to create a seal. The metal-to-metal contact between the cover and body prevents O-ring extrusion.

MAWP (Design Pressure) vs. Temperature Rating for O-ring Models

Material *	Nominal	Order Number	°F	-20 to 100°	150°	200°	250°	300°	350°	400°	450°	500°
(Body, Cover, Seat Ring)	Matorial	Suffix	°C	-29 to 37.8°	66°	93°	121°	149°	177°	204°	232°	260°
Type 316 SS	100 ml 250 ml 500 ml	Add SS 11	psi (bar)	11,000 (758)	11,000 (758)	10,950 (755)	10,900 (752)	10,850 (748)	10,750 (741)	10,650 (734)	10,600 (731)	10,500 (724)
	1000 ml 2000 ml	Add SS 10	psi (bar)	10,500 (724)	10,500 (724)	10,500 (724)	10,450 (721)	10,300 (710)	10,250 (707)	10,050 (693)	10,025 (691)	10,000 (689)
A-286	All Sizes	Add AS 29	psi (bar)	28,800 (1986)	28,500 (1965)	28,450 (1962)	28,200 (1945)	27,850 (1920)	27,350 (1886)	26,500 (1827)	25,700 (1772)	25,200 (1738)
Alloy Steel (SA-723)	All Sizes	Add SL 30	psi (bar)	30,000 (2068)								

Maximum Recommended Operating Pressure (MROP) varies based on gauge, transducer, and rupture disk selection.





^{*} Maximum Allowable Working Pressure

Specifications and Part Numbers


Drawing	Model	Nominal Capacity	Dimensions: Inches (mm)				F375C	F375C		
Number			ID	OD	IL	OAL	Connection Standard	Connection Standard		
30-2557	OR0010SS11 OR0010AS29 OR0010SL30	100 ml	1.25 (32)	2.50 (64)	5.50 (140)	12.75 (324)				
30-2558	OR0025SS11 OR0025AS29 OR0025SL30	250 ml	1.25 (32)	2.50 (64)	13.00 (330)	20.25 (514)	OAL			
30-2528	OR0050SS11 OR0050AS29 OR0050SL30	500ml	2.00 (51)	4.00 (102)	10.87 (276)	17.25 (438)		OAL		
30B-1829	OR0100SS10 OR0100AS29 OR0100SL30	1000 ml	3.00 (76)	6.00 (152)	9.65 (245)	19.80 (503)				
30B-1830	OR0200SS10 OR0200AS29 OR0200SL30	2000 ml	3.00 (76)	6.00 (152)	18.72 (476)	28.80 (732)		OR050/OR100/OR200		

O-ring Closure Compression Cylinders

Models OC and OD O-ring closure compression cylinders are available as single-ended or double-ended vessels. Double-ended vessels allow for easier cleaning. All models have an F375C AE high pressure connection at both ends plus four 3/8"-24 tapped holes for eye bolts and mounting. Standard material of construction is alloy steel. Optional pistons are available. Refer to the model code chart for information on pressure ratings. Other materials of construction and volumes are available on special order.

Specifications and Part Numbers

-								
Drawing		Model	Nominal	Dir				
N	Number	Model	Capacity	ID	OD	IL	OAL	F375C
2	0-9006	0C0100SL15	1000 ml	3.00 (76)	4.94 (125)	9.00 (229)	17.69 (449)	
2	0-9007	0C0100SL30	1000 ml	3.00 (76)	6.00 (152)	9.75 (248)	18.94 (481)	
2	0-9008	0C0100SL60	1000 ml	2.00 (51)	6.00 (152)	19.50 (495)	28.00 (711)	
2	0-9009	0D0200SL15	2000 ml	3.00 (76)	4.94 (125)	18.00 (457)	28.50 (724)	
2	0-9010	0D0200SL30	2000 ml	3.00 (76)	6.00 (152)	18.00 (457)	28.50 (724)	F0750
2	0-9011	0D0200SL60	2000 ml	2.00 (51)	6.00 (152)	39.00 (991)	49.00 (1245)	F375C



Optional O-Ring Chart:

Material	Max. Temperature				
Material	°F	° C			
Nitrile	250°	(121°)			
EPR	300°	(149°)			
Silicone	400°	(204°)			
FKM	450°	(232°)			
PTFE Encapsulated FKM	400°	(204°)			
FFKM	500°	(260°)			

Ordering Guide:



Part Number Example: **OR0025SS11** (example selections indicated in yellow below)

A - Pres	ssure Vessel						
OR	O-Ring Vessel						
ос	O-Ring Compression, Cylinder Single Ended						
OD	O-Ring Closure Compression Cylinder, Double Ended						
B - Cap	pacity						
0010	100 ml						
0025	250 ml						
0050	500 ml						
0100	1000 ml						
0200	2000 ml						
C - Mat	terial of Construction						
SS	316 Stainless Steel						
AS	A-286 Stainless Steel						
SL	Alloy Steel (SA-723)						
D - Pre	ssure: psi (bar)						
10	10,500 (724)						
11	11,000 (758)						
15	15,000 (1734)						
29	28,800 (1986)						
30	30,000 (2069)						
60	60,000 (4137) (OC & OD models only)						

Supporting Information:

Please refer to the following sections of the catalog for complimentary products and additional technical details.

- "MagneDrives®" (Agitation Section, Bulletin AG-SG)
- "Controls" (Instrumentation Section, Bulletin IN-SG)
- "Complete Systems" (Stirred Reactor Section, Bulletin SR-SG)
- "Valves, Fittings, Tubing" (VFT Catalog, www.autoclave.com)

Features and Available Options:

- Different Connections Available
- Other Materials of Construction
- Piston (available on OD series vessels)
- Heaters
- Stands
- ASME Code Stamp, CRN Registration available (consult factory)
- Consult factory to discuss your specific requirements

WARNING

FAILURE, IMPROPER SELECTION OR IMPROPER USE OF THE PRODUCTS AND/OR SYSTEMS DESCRIBED HEREIN OR RELATED ITEMS CAN CAUSE DEATH, PERSONAL INJURY AND PROPERTY DAMAGE.

This document and other information from Parker Hannifin Corporation, its subsidiaries and authorized distributors provide product and/or system options for further investigation by users having technical expertise. It is important that you analyze all aspects of your application and review the information concerning the product or system in the current product catalog. Due to the variety of operating conditions and applications for these products or systems, the user, through its own analysis and testing, is solely responsible for making the final selection of the products and systems and assuring that all performance, safety and warning requirements of the application are met. The products described herein, including without limitation, product features, specifications, designs, availability and pricing, are subject to change by Parker Hannifin Corporation and its subsidiaries at any time without notice.

Offer of Sale

The items described in this document are available for sale by Parker Hannifin Corporation, its subsidiaries or its authorized distributors. Any sale contract entered by Parker will be governed by the provisions stated in Parker's standard terms and conditions of sale (copy available upon request).

©2016 Parker Hannifin Corporation | Autoclave Engineers is a registered trademark of the Parker Hannifin Corporation

06-0055SE

January2016

Autoclave Autoclave Engineers

Instrumentation Products Division

Autoclave Engineers Operation 8325 Hessinger Drive Erie, PA 16509-4679 Tel: 814 860 5700 • Fax: 814 860 5718 www.AutoclaveEngineers.com **Caution!** Do not mix or interchange component parts or tubing with those of other manufacturers. Doing so is unsafe and will void warranty.

Caution! Parker Autoclave Engineers Valves, Fittings, and Tools are not designed to interface with common commercial instrument tubing and are designed to only connect with tubing manufactured to Parker Autoclave Engineers AES specifications. Failure to do so is unsafe and will void warranty.

Bulletin PV-OR

ENGINEERING YOUR SUCCESS.