
ZipperClave®

Stirred Reactors Ordering Guide

500, 1,000, 2,000, and 4,000 ml



Ordering Guide:



| Pressure Vessel | | | | | | |
|-----------------|---|---|---|---|---|---|
| S | S | В | 1 | 2 | 3 | 1 |
| Α | Α | В | С | D | Е | F |

| MagneDrive® Agitator | | | | |
|-------------------------|---|---|---|---|
| Α | 1 | 1 | 1 | A |
| G | Н | J | K | L |

| Internal Accessories | | | | |
|-------------------------|---|---|---|--|
| 1 | 1 | 0 | 1 | |
| М | N | 0 | Р | |

| External Accessories | | | | |
|----------------------|---|---|---|---|
| 1 | D | 1 | 1 | 0 |
| R | S | Т | U | ٧ |

Part Number Example: Z050SSB1231A111A11011D110 (example selections indicated in yellow below)

Model Code

| Volume | |
|--------|------------------------------|
| 050 | 500 ml Zipperclave Reactor |
| 100 | 1,000 ml Zipperclave Reactor |
| 200 | 2,000 ml Zipperclave Reactor |
| 400 | 4,000 ml Zipperclave Reactor |

Pressure Vessel

| AA - Vessel Material | | |
|----------------------|-------------------------------|--|
| SS | 316 Stainless Steel | |
| HC | Hastelloy® C-276 ¹ | |

| B - O-ring Seal Material | | | | |
|--------------------------|--|--|--|--|
| В | Nitrile (Max. Temp. 250°F/121°C) ² | | | |
| С | Ethylene-Propylene (Max. Temp. 300°F/149°C) ² | | | |
| D | PTFE (PTFE Encapsulated FKM) (Max. Temp. 400°F / 204°C) ² | | | |
| E | FKM (Max. Temp. 450°F /232°C) ² | | | |
| F | Silicone (Max. Temp. 400°F /204°C) ² | | | |
| G | FFKM (Kalrez®) 3 (Max. Temp. 500°F /260°C) 2 | | | |

| C - Body Bottom Connection | | | |
|----------------------------|---|--|--|
| 0 | None (No Connection) | | |
| 1 | 1/2" Port Manual valve (requires Floor Stand) 4 | | |
| 2 | AE "Flat Bottom" Connection | | |

| D - Approvals Available 10 | | | |
|----------------------------|-----------------------|--|--|
| 0 | None Required | | |
| 2 | CE Mark and PED | | |
| 3 | Canadian Registration | | |

| E - Stand | | | | |
|-----------|--|--|--|--|
| 0 | None | | | |
| 2 | Tall Bench Top (500 and 1,000 ml ONLY) | | | |
| 3 | Floor | | | |

| F - Body Lift Mechanism | | | |
|-------------------------|------------------|--|--|
| 0 | None | | |
| 1 | Manual Jack | | |
| 2 | Manual Jack (CE) | | |

MagneDrive® Agitator

| | g |
|-------|-------------------------------------|
| G - N | MagneDrive® Agitator |
| Α | MAG075-01 Belt Driven |
| В | iMAG075 Inline |
| С | MAG075-02 Belt Driven |
| Х | No MagneDrive® with opening plugged |

| H - Bearings | |
|--------------|-------------------------------------|
| 0 | None ⁵ |
| 1 | Purebon® 6 (Carbon Graphite) |
| 2 | Fluoropolymer with graphite fiber 7 |
| 3 | Purebon®63310 |

| J - Speed Sensors | |
|-------------------|-----------------------------|
| 0 | None |
| 1 | General Purpose Hall Effect |

| K - N | K - Motors | |
|-------|--|--|
| 0 | None | |
| 1 | DC Variable Speed, 90 VDC, General Purpose | |
| 2 | DC Variable Speed, 180 VDC, General Purpose | |
| 3 | DC Variable Speed, 90 VDC, XP (Non-CE Mark) | |
| 4 | DC Variable Speed, 180 VDC, XP (Non-CE Mark) | |
| 5 | Air with Manual Speed Adjust | |
| 6 | Air with Electronic Speed Adjust | |
| 7 | AC Motor, XP CE Mark | |
| С | Belt & Guard WITHOUT MOTOR | |
| D | 1/8 HP 0-130 VDC Variable Speed GP Inline | |
| E | 1/3 HP 0-130 VDC Variable Speed GP Inline | |
| F | Air Motor - Manual Speed Adjust Inline | |
| G | Air Motor - Electronic Speed Adjust Inline | |

| L - Impellers / Shaft / Baffles | | |
|---------------------------------|--|--|
| Α | AE Dispersimax™ (6 blades) with Baffle Bar | |
| В | Turbine (6 blades) with Baffle Bar | |
| С | Axial-Up (4 blades) with Baffle Bar | |
| D | Axial-Down (4 blades) with Baffle Bar | |
| Х | None ⁵ | |





Internal Accessories

| M - Liq | M - Liquid Sample | |
|---------|--|--|
| 0 | None, Plugged Connection | |
| 1 | Sample Tube Only | |
| 2 | Sample Tube with Manual Valve | |
| 5 | Sample Tube with Manual Valve & Filter | |

| N - Blow Pipe | |
|---------------|-----------------------------|
| 0 | None, Plugged Connection |
| 1 | Blow Pipe Only |
| 2 | Blow Pipe with Manual Valve |

| - | O - Sparge Tube | |
|---|-----------------|-------------------------------|
| ľ | O - Sparge Tube | |
| | 0 | None, Plugged Connection |
| | 1 | Sparge Tube Only |
| Г | 2 | Sparge Tube with Manual Valve |

| P - Cooling Coil | |
|------------------|---|
| 0 | None, Plugged Connection |
| 1 | Cooling Coil Only |
| 2 | Cooling Coil with Manual Valve |
| 3 | Cooling Coil with Solenoid Valve (120 Volt) |
| 4 | Cooling Coil with Solenoid Valve (240 Volt) |

- 1. HASTELLOY® is a registered trademark of Haynes International Inc
- 2. Temperature limits are suggested. Actual performance will vary with chemical
- 3. Kalrez® is a registered trademark of DuPont.
- 4. The drain valve is a "Flush" design (no dead volume) that extends approximately 8.25" (210 mm) below the vesssel.
- 5. Use this option only if X (No MagneDrive®) is selected as the model of MagneDrive® agitator
- 6. Purebon® is a registered trademark of Morgan AM&T.
- 7. Fluoropolymer bearings have a maximum recommended service temperature of 500°F (260°C).
- 8. MROP may be further reduced by temperature and number of cycles.
- 9. When heating/cooling is selected, the reactor is supplied with a process Type K Thermocouple and Thermowell, and an external Type K Thermocouple. When no heating/cooling is selected, the reactor will be supplied with a plugged connection for the process thermocouple.
- 10. Consult factory for pricing and rating of code vessels

External Accessories

| R - Ven | R - Vent Valve | |
|---------|--|--|
| 0 | None, Plugged Connection | |
| 1 | Vent with Manual Valve | |
| 2 | High Volume Vent with Solenoid Valve (120 Volt) | |
| 3 | High Volume Vent with Solenoid Valve (240 Volt) | |
| 4 | BPR Digital (120 Volt) | |
| 5 | BPR Digital (240 Volt) | |
| 7 | BPR Digital with High Volume Vent 120 VAC Solenoid | |
| 8 | BPR Digital with High Volume Vent 240 VAC Solenoid | |

| S - Pressure Gauge/Transducer (MROP = Max. Recommended Operating Pressure) | |
|--|---|
| Α | 600 psi Gauge Only (450 psi MROP) 8 |
| В | 1,000 psi Gauge Only (750 psi MROP) ⁸ |
| С | 2,000 psi Gauge Only (1,500 psi MROP) 8 |
| D | 3,000 psi Gauge Only (1,880 psi MROP) 8 |
| G | 600 psi Gauge/1 ksi Transducer (450 psi MROP) 8 |
| Н | 1,000 psi Gauge/1 ksi Transducer (750 psi MROP) 8 |
| J | 2,000 psi Gauge/3 ksi Transducer (1,500 psi MROP) 8 |
| K | 3,000 psi Gauge/3 ksi Transducer (1,880 psi MROP) 8 |
| N | 600 psi Gauge/1 ksi IS Transducer (450 psi MROP) 8 |
| Р | 1,000 psi Gauge/1 ksi IS Transducer (750 psi MROP) 8 |
| Q | 2,000 psi Gauge/ 3 ksi IS Transducer (1,500 psi MROP) 8 |
| R | 3,000 psi Gauge/3 ksi IS Transducer (1,880 psi MROP) 8 |

| T - Heating and Cooling 9 | | |
|---------------------------|---|--|
| 0 | None | |
| 1 | Electric 120 VAC, Single Phase | |
| 2 | Electric 240 VAC, Single Phase | |
| 3 | 120 VAC, Purgeable Furnace | |
| 4 | 240 VAC, Purgable Furnace | |
| 5 | Baffled Removable Jacket, 1/4" FNPT Connections 450°F (232°C) Maximum | |

| U - Gas Inlet | |
|---------------|--|
| 0 | None, Plugged Connection |
| 1 | Gas Inlet Line with One (1) Manual Valve |
| 2 | Gas Inlet Line with Two (2) Manual Valve (Shared Connection) |
| 3 | Forward Pressure Regulation (FPR) - Digital 120VAC |
| 4 | Forward Pressure Regulation (FPR) - Digital 240 VAC |

| V - Charging Valve | |
|--------------------|---|
| 0 | None, Plugged Connection |
| 1 | 3/8" Manual Charging Valve |
| 2 | Manual Valve with 8cc Charging Cartridge |
| 3 | Manual Valve with 20cc Charging Cartridge |
| 4 | Reflux Condenser |

WARNING

FAILURE, IMPROPER SELECTION OR IMPROPER USE OF THE PRODUCTS AND/OR SYSTEMS DESCRIBED HEREIN OR RELATED ITEMS CAN CAUSE DEATH, PERSONAL INJURY AND PROPERTY DAMAGE. This document and other information from Parker Hannifin Corporation, its subsidiaries and authorized distributors provide product and/or system options for further investigation by users having technical expertise. It is important that you analyze all aspects of your application and review the information concerning the product or system in the current product catalog. Due to the variety of operating conditions and applications for these products or systems, the user, through its own analysis and testing, is solely responsible for making the final selection of the products and systems and assuring that all performance, safety and warning requirements of the application are met. The products described herein, including without limitation, product features, specifications, designs, availability and pricing, are subject to change by Parker Hannifin Corporation and its subsidiaries at any time without notice.

Offer of Sale

The items described in this document are available for sale by Parker Hannifin Corporation, its subsidiaries or its authorized distributors. Any sale contract entered by Parker will be governed by the provisions stated in Parker's standard terms and conditions of sale (copy available upon request).

©2015 Parker Hannifin Corporation | Autoclave Engineers is a registered trademark of the Parker Hannifin Corporation



November2015





Instrumentation Products Division

Autoclave Engineers Operation 8325 Hessinger Drive Erie, PA 16509-4679 Tel: 814 860 5700 • Fax: 814 860 5718 www.AutoclaveEngineers.com

Caution! Do not mix or interchange component parts or tubing with those of other manufacturers. Doing so is unsafe and will void warranty.

Caution! Parker Autoclave Engineers Valves, Fittings, Systems, and Tools are not designed to interface with common commercial instrument tubing and are designed to only connect with tubing manufactured to Parker Autoclave Engineers AES specifications. Failure to do so is unsafe and will void warranty.

